


Carl Bermingham

# Lecture 6: Cleaning & Sterilisation Validation






# MCQ 2

- The second MCQ for your Biopharmaceutical Validation module will take place on **12<sup>th</sup> March 2025**.
- There will be **20 multiple choice questions** (one correct answer for each) and once you begin your attempt you will have **50 minutes to complete**.
- The questions will cover **lectures 4-6**. The MCQ is worth **10% of your overall grade**.
- You will have **one attempt** only and the 50 minutes will begin when you click "start attempt".
- The MCQ will be on Moodle. **It will be open from 7.30am on 12<sup>th</sup> March to 7.30am on 13<sup>th</sup> March 2025**.
- **IMPORTANT:** attempts outside of this time frame will not be accepted except in certain extenuating circumstances. This is to ensure fairness to all students.



# Topics





# Cleaning in the Biopharma Industry

## Why is it necessary?

- Biopharmaceuticals are administered intravenously
- Risk of contamination must be reduced
- Sterility & cleanability become the highest priority

## Regular Cleaning of Equipment

- Reduces and controls bioburden
- Prevents cross contamination
- Sets baseline for subsequent steaming and sterilising


# Types of Cleaning



<https://amcrasto.wordpress.com/tag/cleaning-validation/>

1. Manual Cleaning:
  - variation between operators

2. Automated Cleaning
  - can be validated  
i.e. less variation






# Cleaning Strategies

- **Clean in Place (CIP)**

- Equipment stays in place (at point of process use) and cleaning solutions are routed through hoses and piping from Utilities to the equipment.
- E.g. bioreactors and large immobile vessels



- **Clean out of Place (COP)**


- Partswasher: Similar to a dishwasher used in homes. Parts and glassware are loaded into a chamber and solutions are sprayed onto the parts.
- COP stations for mobile vessels
- Sonicating baths e.g. some sparger heads








# CIP of Vessel





# Mains/Pipework CIP






## Design considerations:


- Split routes – total route coverage
- Avoid Dead ends
- Avoid tees
- Refer to ASME BPE guide

## Standard is 1.5-2.1 m/s velocity

- Important to ensure turbulent flow for flat velocity profile and thin boundary layer
- Turbulent flow minimises ineffective movement and cleaning in boundary layer.
- Turbulent flow also cleans out dead legs, sample points and instrumentation better than laminar profile – These elements are still preferably avoided



If Deadlegs are unavoidable -  
Ensure flow directed into dead end  
and is as short as possible - L/D = 1.5







# Cleaning Validation

- There are two main factors which affect the cleaning of a vessel:
  1. The physical design of both the vessel and the associated cleaning system/components.
  2. The cleaning cycle and its parameters
- Both of these factors need to be considered and validated appropriately .



# Physical Design: Spray Ball Devices

## Spray device selection:

- deliver solutions to all locations
- vessel size
- vessel geometry
- agitator, spargers, dip tubes, sample ports
- shadowing

## Spray ball can be:


- fixed or rotating





# Physical Design: Coverage Testing

- Carried out during **FAT/OQ**
- Tests spray ball functionality
  - i.e. ability to deliver solution to all parts of vessel
  - Coverage test is not a demonstration of how clean the vessel is, but how cleanable it is with the current spray ball design
- Fluorescent solution (riboflavin) is sprayed on all internal surfaces.
- A water cycle is run through the spray ball
- Once the cycle is complete, the vessel is opened, and all internal surfaces are inspected using a UV light
- If the spray ball design and location is effective, all riboflavin should have been removed during the water cycle, and not fluorescence should be observed.



Internal surfaces of a vessel covered with Riboflavin

## What conditions?

- short (5 mins), ambient water rinse



# Coverage Test Failure – e.g. FAT, OQ

- Incorrect design of the spray ball
- Incorrect design of lines or ports of the vessel
- Grease or other residues in the vessel
- Spray ball not installed correctly
- Supply water pressure not sufficient
- Blockage of spray ball

Identify & fix root cause – repeat test



# Cleaning Cycle Development

- Different recipes tested
- Effective cleaning cycle is developed which
  - Removes all residues
  - Leaves equipment dry
- Least amount of water
- Coolest temperature
- Shortest amount of time



Saves money





# Cleaning Cycle Parameters

Ideally, cleaning cycles should optimise the parameters below, to generate the most efficient cleaning cycle possible in terms of cost, energy and efficacy.



## TEMPERATURE

Some chemicals and flushes work better at hot temps, but can a cycle be completed as effectively at 60°C rather than 80°C?



## CHEMICAL CONCENTRATION

Need to use the optimal concentration of chemical to effectively clean the surface, but is it preferable and safer to use a chemical at the lowest effective concentration?



## FLOW RATE & PRESSURE

High flow rates and pressurised liquid streams can physically remove soils, but ideally, we want to use as little water as possible – time and energy!



## CONTACT TIME

We need to allow chemicals time to “work”, but ideally, this should be as short a time as possible. Longer contact time is necessary at lower temps and lower concentration.



# Cleaning Validation

- Must be performed *“to confirm the effectiveness of any cleaning procedure for all product contact equipment”*  
EUGMP Vol 4. Annex 15
- Once the cleaning cycle has been developed, it needs to be formally validated.
- Can be performed before, or in tandem with Process Validation (PV).
- Consumes a considerable amount of time and resources.
- Validation testing to demonstrate cycle efficacy include:

## Visual Inspection

- Riboflavin Testing – coverage test.
- Visual “dryness” check after cycle completion.

## Indirect Sampling

- also known as rinse sampling, performed on the WFI rinse prior to blow-down. Samples taken for TOC, conductivity, bioburden and endotoxin.

## Online testing

- Examples of PAT – online conductivity and TOC meters

## Direct sampling

- Swab sampling for TOC and bioburden





# Validation - Acceptance Criteria for Dryness

- Dry, no water droplets visible
- Isolated water droplets visible
- Numerous water droplets visible, no pooling
- Pooling evident at bottom of vessel



<https://ispe.org/pharmaceutical-engineering/september-october-2017/proof-closure-life-cycle-closed-systems>



# What to Sample For?

1. Conductivity

2. Microbes - Endotoxin

3. Protein - Drug substance

4. TOC – Total Organic Carbon





# Clean/Dirty Hold Times

## NIBRT:

Clean hold time:  $\leq 14$  days

Dirty hold time:  $\leq 3$  days

Steam hold time:  $\leq 24$ hrs

## Vessel Cleaned

Vessel Name: \_\_\_\_\_

Vessel No.: \_\_\_\_\_

COP Date: \_\_\_\_\_

Expiry of COP State: \_\_\_\_\_


Prepared By: \_\_\_\_\_

Date: \_\_\_\_\_

- How are these validated?
- What happens if hold times are exceeded?



# Topics





# Sterilization vs. Depyrogenation

- Sterilization: Validated process used to render a product free of living microorganisms including bacterial endospores.
- Depyrogenation: Removal or inactivation of bacterial endotoxin. Endotoxins are toxins present within certain bacterial cell walls, which are released when the cell is broken down.
- An overkill approach is applied when sterilizing - assurance of  $10^6$  probability of non-sterility.



# Sterilisation Methods

- **Steam** -  $\sim 121^{\circ}\text{C}$  – shorter exposure times, higher kill rate, more complicated system design
- **Dry Heat** -  $\sim 400^{\circ}\text{C}$  – longer exposure times, lower kill rate, simpler system design
- Vaporised Hydrogen Peroxide (VHP) – ok for sterilizing background work surfaces (e.g. internal surfaces of filling line) but not suitable for direct and indirect product contact parts.
- Radiation - typically used for single use plastics



# 1. Dry Heat Sterilization & Depyrogenation

- Dry heat is the method of choice for sterilizing items which will tolerate high temperatures
- Generally less complicated than steam processes
- Higher temperature and/or longer exposure times are required
- Microbial lethality associated with dry heat is much lower than that for saturated steam at the same temperature





# Depyrogenation

- Refers to the removal of pyrogens – defined as “any substance that can cause a fever. Bacterial pyrogens include endotoxins and exotoxins”
- There is no uniform regulation in relation to the time and temperature required for depyrogenation
- 3 log step reduction of endotoxins is required by the USP and FDA
- No European regulation for depyrogenation of final containers for parenterals



## 2. Steam Sterilization - Steam In Place (SIP)


- Widely adopted method of in-line sterilisation for process equipment.
- Very effective if done right – needs to be considered at early concept stage and considerations incorporated into facility design.
- Critical design parameters:
  - Proper steam distribution and contact with all surfaces
  - Non-condensable gas removal – Steam traps.
  - Continuous condensate removal.
- Also need to consider vent filters for gas inlets/outlets – must be installed before steaming – therefore must be suitable for steaming – must remain dry for batch
- Use hydrophobic materials – prevent blockage by humidity



# Steam In Place

## Moist heat sterilisation


- Clean steam
- Uses latent heat of steam
- Contacts a surface and gas condenses releasing heat to the surface
- Microbes destroyed on the surface





# SIP Consideration


- Large amounts of Steam are needed per Cycle
  - Steam vessel first – Vent filter limits flowrate
- Appropriate placement and usage of steam traps
  - Remove non-condensable gases
- Correctly sloped/placed drains
  - Ensure adequate condensate removal
- Efficient SIP procedure and cycle
  - Test pressure, validation study and post SIP assurance






# Air Removal

- Air is a mixture of non-condensable gas
- Air is an excellent insulator
- A small layer can prevent steam accessing a surface
- ∴ As much air must be removed as possible






# What you don't want: Bi-Directional Steam





# What you do want: Uni-Directional Steam



Saturated steam introduced from one direction.  
Air gets pushed to multiple outlets





# SIP Cycle

## Pre-Cycle

- Pressure hold test,
- Flush out non-condensable gases and heat system


## Steam Cycle

- Continues for time based on validation study
- Temperature and pressure are monitored to maintain saturated steam conditions
- Condensate continuously drains due to system design and is replaced with fresh steam
- Air is removed continuously through bleed valves and steam traps

## Post SIP


- Venting with sterile gas (typically process air) to remove condensate
- Retain positive pressure during cooling to avoid vacuum conditions
- May choose to integrity test vent filter for assurance

# Sterilisation Plateau (121°C)






# Filter Drying and Cooldown





# Awaiting Loading





# Steam Sterilization Parameters & Principles






# Steam Sterilization Validation

## 1. Thermocouples

- Verify correct temperature, pressure, and time are reached



## 2. Biological Indicators

- Verify all organisms are killed



# Thermocouples

- Two dissimilar metals – insulated wires
- Hot junction / cold junction
- Generates EMF within the wires
- This current is dependent on temperature
- Current EU and US regulatory requirements - EN285, HTM2010, HTM 01-05, EN17665, PDA Technical Report 01.







# Biological Indicators

- **Contain bacterial endospores**
  - *Geobacillus stearothermophilus*
  - *Bacillus atrophaeus*
- **Spore strips**
  - inoculated with  $10^6$  spores
- **Incubated after steaming cycle**
- **Some BIs contain growth medium & colour indicator**





# Which Locations?

- Worst case scenarios
- Coldest point? – steam rises, condensate drains
- Areas of difficult steam penetration
- Places that may trap air pockets

### 3. Autoclave Steam Sterilization – Steam Out of Place (SOP)



- Smaller equipment – hard/wrapped items or liquid load
- Shorter times and lower temperatures than dry sterilization
- Also uses latent heat of steam – direct steam contact is required
- Cost-effective
- A dedicated pure steam supply is recommended



# Autoclave Cycle

- **Pre-Conditioning:** Air is removed, and load is humidified and heated
- **Exposure:** Temperature is raised, and Exposure time is conducted
- **Post-Conditioning:** load is cooled and dried, chamber is brought to atmospheric pressure
- Cycles will vary with individual applications – vital to ensure a suitable cycle is being used



# Autoclave Validation Tests

- Similar to the SIP of vessels and larger systems, autoclaves will be tested during validation studies (and periodic re-validation studies) using thermocouples and biological indicators.

The following tests are also typically carried out on a regular basis post-validation as part of routine monitoring:

- Leak rate test (Before use) – check that autoclave is sealed and has no leaks
- Bowie-Dick (Daily) – check penetration power of steam
- Air detector function (Weekly) – ensure air detector functionality
- Monitoring of chamber and printouts



- Routine monitoring tests must be validated as fit-for-use/suitable during system development/validation.



# Thank You

