

Level 8 Cell Culture Processing (BIO08045)

Lecture 3 – “Cell Banking & Cell Culture Media Design”

- Dermot O’ Sullivan
- Dermot.osullivan@nibrt.ie

Learning Objectives

The importance of cell banks

Cell line verification

Describe the key components in animal cell culture media

Compare and contrast serum-supplemented, serum-free media and protein free media formulations


Recommended Reading

- R.I. Freshney, *“Culture of Animal cells : a manual of basic technique”*. 8th edition, Wiley-Blackwell (2021)
- John M. Davis, *“Animal Cell Culture: Essential Methods”*, 1st Edition. Wiley, (2011)
- Michael Butler, *“Animal Cell Technology”*. 2nd edition, BIOS Scientific (2004)
- JM Davis, *“Basic Cell Culture: A Practical Approach”*. 2nd edition, Oxford University Press (2002)
- Terence Cartwright, *“Animal Cells as Bioreactors”*, Cambridge University Press, 1994

Websites/Videos

- Some useful videos of cells in culture can be viewed at:
<https://www.sartorius.com/en/products/live-cell-imaging-analysis>
- American tissue culture collection : www.atcc.com
- Fundamental techniques in Cell Culture (a laboratory handbook): available from Sigma website (in pdf or hard copy)
- Cells alive website (<http://www.cellsalive.com>) particularly for those with little biology

Lecture Topics



Cell Line & Process Development

Early stage



- Cell line development
- Production media selection
- Clone with the highest productivity properties selection

Late stage




- Feeding and supplementation strategy optimization
- Process parameters optimization
- Product titer and quality optimization

Process transfer



- Three consolidation runs in R&D before transfer
- Data alignment between R&D and Production Units


Cell Line & Process Development



- Knowing gene for the protein you want is great, but what cell line to use? What clone form that cell line is best. 100s of possibilities!
- 60 or more nutritional components in culture media, how many combinations? When to feed them? Inducers, promoters?
- What temperature? What oxygen level? CO₂? pH any shifts? When to harvest?
- A strategy of multi-factorial design is the natural way to attack this type of problem, but is difficult to execute in cell culture because the parameters interact strongly-requiring a lot of experiments. This means models!

Cell Line Development to Cell Banking

From: Comprehensive characterization of glutamine synthetase-mediated selection for the establishment of recombinant CHO cells producing monoclonal antibodies Soo Min Noh¹, Seunghyeon Shin¹ & Gyun Min Lee.
Scientific REPOrTS | (2018) 8:5361 | DOI:10.1038/s41598-018-23720-9



A schematic diagram of the process for mAb producing clone generation and long-term culture for testing the production stability.

<https://www.nature.com/articles/s41598-018-23720-9/figures/1>

Cell Banking

- Secure, controlled long term storage of cell stocks:
 - Required to control cell lines, to have original stocks for each batch or in case of contamination.
- Required for regulation and validation of process:
 - Common starting point for each batch of product.
- Frozen stocks of cells reduces risk of contamination:
 - Stocks are identical avoiding genetic instability in cell lines.
- Vials stored in a cryofreezer or in liquid nitrogen:
 - Cells frozen in 1-2mL cryovials.
 - Sub zero temperatures are used which effectively stop any biological activity that could lead to cell death.
 - A cryoprotective agent is used to prevent damage to cells by the formation of crystals on freezing. A protective agent like DMSO (Dimethyl sulphoxide) is added to the cells before freezing.




Cell Banking

- Cell banks are mostly stored in liquid nitrogen tanks/dewars:
 - Cells can be stored in the liquid or in the vapour of the liquid nitrogen tank.
 - Temperature of liquid nitrogen -196°C .
 - ✓ Constant ultra low temperature.
 - ✗ Contamination is possible if vials are not secured.
 - Temperature of liquid nitrogen vapour -130°C .
 - ✓ Less risk of cross contamination.
 - ✗ Poor temperature gradient.



Cell Banking



- Ultra-low temperatures.
- Stops biological activity that could lead to cell death.

Cell Banking



- Stocks prepared from primary culture-early passage numbers:
 - Cell bank has a 2 tier system which consists of **master cell bank (MCB)** and **working cell bank (WCB)**.
 - Master stocks are prepared for the **master cell bank**.
 - From this, stocks are prepared for a **working cell bank**.
- Vials are thawed for consistency checks:
 - Cell banks can be stored in several locations.



Remove Vial from the Working Cell Bank and expand culture for each production batch



Cell Banks



Principles of Good Cell Culture Practice


- In all aspects of sourcing, banking and preparing cell cultures, the principles of Good Cell Culture Practices (GCCPs) should be observed. Fundamental features to be considered in the development of cell cultures for production or testing are:

1. **Authenticity**, including identity, provenance, and genotypic/phenotypic characteristics;
2. **Absence of contamination** with another cell line;
3. **Free of microbiological contamination**, and
4. **Stability** and functional integrity on extended in vitro passage.

Cell Line Characterisation

- What is really in the vessel? MCB and WCB are among the major assets of a company.
- Need to ensure that you have good quality: check number of parameters:
 1. **Viability**: when you set up MCB must thaw vial to check viability;
 2. **Contamination**: checks for bacteria, fungi and especially for mycoplasma;
 3. **Morphology**: shape cells, attachment;
 4. **Isoenzymes**: can be used to check for species e.g. Lactate dehydrogenase;
 5. **Karyotype**: chromosomal analysis;
 6. **DNA fingerprinting**: more expensive and usually done by specialised laboratory;
 7. **Productivity/Titre**: check levels of product.



Cell Line Characterisation




Light microscope

Karyotype

*Spectral Ray
Computer image*



DNA fingerprint



Isoenzyme analysis

Cell Line Characterization Studies

	MASTER CELL BANK	MANUFACTURERS WORKING CELL BANK	END OF PRODUCTION CELLS <small>Cells at the limit of in vitro cell age</small>
Authenticity & Characteristics			
Identity	Yes	Optional	Yes
Growth / Morphology	Yes	Yes	Yes
Genetic Stability	Program Dependent	Program Dependent	Program Dependent
Tumorigenicity	Cell Line Dependent	No	No
Adventitious Agents			
Sterility	Yes	Yes	Yes
Mycoplasma	Yes	Yes	Yes
General Virus Screens	Yes	No	Yes
Species-Specific Virus Tests	Yes	No	No
Retrovirus Tests	Yes	No	Yes
Adventitious Agents from Raw Material	Yes	Program Dependent	Program Dependent


The table above provides an overview of the tests typically required to characterize mammalian cell lines as outlined in the worldwide regulatory guidance documents - which tests to conduct and at which stage.

Cell Line Characterization

- See Recommendations for the evaluation of animal cell cultures as substrates for manufacture of biological medicinal products and for the characterization of cell banks. *WHO*, 2010.

<https://www.who.int/publications/m/item/evaluation-of-cell-substrates-for-the-production-of-biologicals-revision-of-who-recommendations>

Lecture Topics




Cell Culture Media

- Cell culture media must contain all the components necessary to support the growth of cells.
- Components must be present in the correct form and concentration.

Carbon/energy source + nitrogen
source + other nutrients

=

Cell biomass + product + CO₂ + H₂O +
heat + metabolites



History of Cell Culture

- Pioneer scientists such as Dulbecco, Ham & Eagle determined the nutrient requirements of animal and human cells in the 1950s.
 - They did so by measuring nutrient uptake of these cell lines, and also by trial & error.
- **DMEM-F12** is one of the most popular base media, it is a 50/50 mixture of Dulbecco's Modified Eagle's Medium & Ham's Nutrient Mixture F12 (see following table for the constituents).

Cell Culture Media

- Medium optimization is a large part of the bioprocess and depends on many factors:
 1. Feeding pattern-batch, fed batch, perfusion
 2. Bioreactor design
 3. Culture conditions-cell type, density
 4. Preparation format-powder, liquid
 5. Downstream purification
 6. Product application
 7. Business factors including cost, equipment, yield etc

Chemical environment of cells

- Growth media contains 2 parts:
 1. Complex / Non-defined portion : usually serum
 2. Defined part : glucose, amino acids, vitamins, salts etc
- First media developed by Eagle – called Eagle’s minimal essential medium (**EMEM**)
- Later modified by Dulbecco – referred to as **DMEM** (Dulbecco’s modification of Eagle’s MEM)
- Other common media such as **HAMs F-12, M199**

*For lab studies.
In biomanufacturing today,
predominantly serum-free,
protein-free*

Cell Culture Media Components

There are two main types of media used in cell culture processes:

Complex Media:		Defined Media:	
Broad range of nutrients, vitamins, minerals etc.	Adventitious agents and other contaminants.	Defined concentration of components	Potentially lower growth
Culture of many different cell types	Batch to batch variation	Little batch to batch variation	Cells may need to be adapted to grow without serum
High cell density and productivity	Supply may be restricted due to disease or drought	Consistent quality	Optimisation of formula can be difficult
Shear protection	Difficulties in downstream processing due to large quantities of proteins etc.	Lower risk of contamination	Higher likelihood of shear stress
Less expensive	Disliked by regulators	Preferred by regulators	Expensive

Lecture Topics

A horizontal bar with a grey circle on the left. A line connects the top of the circle to the top-left corner of the bar.

Cell Banking

A horizontal bar with a grey circle on the left. A line connects the top of the circle to the top-left corner of the bar.

Cell Culture Media

A horizontal bar with a blue circle on the left. A line connects the top of the circle to the top-left corner of the bar.

Defined Portion of Cell Culture Media

A horizontal bar with a white circle on the left. A line connects the top of the circle to the top-left corner of the bar.


Complex Portion of Cell Culture Media

Defined Portion of Media

Carbon/ Energy source	<ul style="list-style-type: none"> • Carbohydrate (usually glucose)
Nitrogen source	<ul style="list-style-type: none"> • Ammonia and amino acids (Amino acid and protein synthesis)
Buffers	<ul style="list-style-type: none"> • For pH control (often bicarbonate buffer)
Vitamins	<ul style="list-style-type: none"> • Riboflavin, thiamine, ascorbic acid etc.
Growth Factors	<ul style="list-style-type: none"> • Signalling hormones (IGF-1, IL-6, G-CSF etc.)
Inorganic salts/ Trace metals	<ul style="list-style-type: none"> • Sodium chloride (osmolality) • Iron, zinc, magnesium, potassium, calcium etc.

Media pH

- Need to maintain a pH typically between 6.8 - 7.5 for mammalian cell cultures
 - Use phenol red as an indicator; media is initially red and changes to yellow as cells grow and metabolise
- Cytochromes – require neutral pH
 - Essential in cellular energy (ATP) production



DMEM as an example

Media formulations		DMEM	DMEM with stable glutamine
Components		130-091-437 g/L	130-091-438 g/L
→	Inorganic salts		
	CaCl ₂ · 2H ₂ O	0.264	0.264
	Fe(NO ₃) ₃ · 9H ₂ O	0.001	0.001
	KCl	0.4	0.4
	MgSO ₄ · 7H ₂ O	0.2	0.2
	NaCl	6.4	6.4
	NaHCO ₃	3.7	3.7
	NaH ₂ PO ₄ · H ₂ O	0.125	0.125
→	Other components		
	Glucose	4.5	4.5
	Phenol red Na	0.0159	0.0159
→	Sodium Pyruvate	0.11	0.11
	Amino acids		
→	L-Arginine.HCl	0.084	0.084
	L-Cystine	0.048	0.048
	L-Alanyl-L-Glutamine	—	0.868
	Glycine	0.03	0.03
	L-Histidine HClH ₂ O	0.042	0.042
	L-Isoleucine	0.105	0.105
	L-Leucine	0.105	0.105
	L-Lysine.HCL	0.146	0.146
	L-Methionine	0.030	0.030
	L-Phenylalanine	0.066	0.066
	L-Serine	0.042	0.042
	L-Threonine	0.095	0.095
	L-Tryptophan	0.016	0.016
	L-Tyrosine	0.072	0.072
	L-Valine	0.094	0.094
	Vitamins		
	D-Ca Pantothenate	0.004	0.004
	Choline Chloride	0.004	0.004
	Folic acid	0.004	0.004
	myo-Inositol	0.0072	0.0072
	Nicotinamide	0.004	0.004
	Pyridoxine.HCl	0.004	0.004
	Riboflavin	0.0004	0.0004
	Thiamine.HCl	0.004	0.004
Reference:			
1. Dulbecco et al. (1959) Virology 8: 396-397.			

DMEM-F12 Components: Inorganic Salts

Component Amount (mg/L)

• CaCl_2	116.6
• CuSO_4	0.001
• $\text{Fe}(\text{NO}_3)_3 \cdot 9\text{H}_2\text{O}$	0.05
• $\text{FeSO}_4 \cdot 7\text{H}_2\text{O}$	0.42
• MgCl_2	28.6
• MgSO_4	48.4
• KCl	312
• NaHCO_3	2438
• NaCl	7000
• Na_2HPO_4	71
• NaH_2PO_4	62.5
• $\text{ZnSO}_4 \cdot 7\text{H}_2\text{O}$	0.43

- **Provide metals needed for cell metabolism:**
 - Included as salts (e.g. chlorides & sulphates for improved solubility).
- **DMEM-F12 is low in iron.** Therefore, supplement with extra iron (Ferric Ammonium Sulphate or Ferrous Citrate, up to 100 mg/L).
- Sodium bicarbonate (NaHCO_3) - buffering system to maintain pH. In bioreactors use CO_2 sparging to complete the buffer system.
- Sodium Chloride (NaCl) conc. is often lowered to compensate for the extra supplements added to DMEM-F12 in order to maintain an osmolarity of around 320mOsm/Kg

DMEM-F12 Vitamins

- **Vitamin Amount (mg/L)**
 - Biotin 0.003
 - Ca Pantothenate 2.2
 - Choline.Cl 9
 - Cyanocobalumin 0.68
 - Folic Acid 2.65
 - I-Inositol 12.6
 - Niacinamide 2
 - Pyridoxine.HCl 2
 - Riboflavin 0.2
 - Thiamine.HCl 2.2
- Vitamins are included because many enzymes need co-factors (such as vitamins and metal ions) to catalyze biochemical reactions.

DMEM-F12 Amino Acids

- All 20 common amino acids in L-form (mg/L):


○ Alanine	4.5	○ Isoleucine	54.5
○ Arginine.HCl	147.5	○ Leucine	59
○ Asparagine.H ₂ O	7.5	○ Lysine.HCl	91.3
○ Aspartic Acid	6.7	○ Methionine	17.2
○ Cysteine.HCl.H ₂ O	17.6	○ Phenylalanine	35.5
○ Cystine.2HCl	81.3	○ Proline	17.3
○ Glutamic Acid	7.4	○ Serine	26.3
○ Glutamine	365	○ Threonine	53.5
○ Glycine	18.8	○ Tryptophan	9
○ Histidine.HCl.H ₂ O	31.5	○ Tyrosine.Na ₂ .2H ₂ O	55.8
		○ Valine	52.9

DMEM-F12 Amino Acids

- **Arginine and Glutamine** are used up quickly by the cell because they feed directly into the Tri-Carboxylic Acid cycle (TCA, or Krebs cycle) and are an alternative source of energy, in addition to their structural role.
 - Many amino acids are inter-convertible through transaminase reactions.
- Glutamine is unstable in solution & breaks down to generate ammonia (an unwanted metabolite).
 - Therefore glutamine (along with glucose) is often included in bioreactor feed regimes.

Interplay of Nutrients in CHO Cell Status


From: Dickson, A.J. Enhancement of production of protein biopharmaceuticals by mammalian cell cultures: the metabolomics perspective. *Current Opinion in Biotechnology* 2014, 30:73–79 <http://dx.doi.org/10.1016/j.copbio.2014.06.004>



Key metabolic indicators of CHO cell function. This figure summarises the metabolites that have been defined from a variety of studies as potential indicators of CHO cell function (in terms of growth and recombinant protein production). Glucose, glutamine and asparagine are shown as the major carbon sources utilised during cell culture and other metabolites shown in black text present metabolites that accumulate during cell culture (with different timing or phases of appearance associated with the status of cell processes). Metabolites in red are produced in limited amounts at late stages (stationary and onwards) but provide markers of CHO cell status and phenotype (in both intracellular and or extracellular localisation). The intracellular metabolites shown in red present a consensus range of metabolites associated with the quality of CHO status during bioprocesses (as described within this review).

Cell Metabolism

http://upload.wikimedia.org/wikipedia/en/b/b4/Citric_acid_cycle_noi.GIF



DMEM-F12 Other Components

- Other Components (mg/L)

○ Glucose	2151	} Energy
○ Na.Pyruvate	55	
○ Na.Hypoxanthine	2.4	} Precursors for DNA synthesis
○ Thymidine	0.37	
○ Linoleic Acid	0.04	} Essential fatty acids
○ Lipoic Acid	0.1	
○ Putrescine.2HCl	0.08	
○ (Phenol Red)	(8.6) – used in laboratory media as a pH indicator	

DMEM-F12 Other Components

- **Glucose is the main energy source for the cell.**
 - However, CHO cells are not good at completing glucose metabolism to generate the 38 ATP molecules of full oxidative metabolism.
 - They often stop at glycolysis, yielding lactate.
 - Hence amino acids are supplemented to provide other energy sources.
- Sodium pyruvate readily enters the TCA cycle.
- Arginine and Glutamine feed directly into the Tri-Carboxylic Acid cycle (TCA, or Krebs cycle) and are an alternative source of energy

Other Supplements & Notes (1)

- Base medium is often supplemented with insulin (10mg/L) to increase glucose uptake and promote cell growth (putrescine does this also).
- Plant or yeast hydrolysates (extracts that are digested by proteases or acid hydrolysis) are a good source of small peptides, which are often taken up by cells more easily than free amino acids. These are added at 5-10 g/L.
- **Trace elements** such as vanadium, cobalt & iodide are often supplemented in very small amounts (<0.0001 mg/L) to DMEM-F12 as they are co-factors to some enzymes. However, they are toxic to cells at higher doses

Other Supplements & Notes (2)

- Today, most media used in manufacturing do not contain animal-derived components such as serum, albumin or transferrin, for safety reasons.
 - Chemically defined and protein free media – more later!
 - In the lab, media often come as sterile liquids
- At manufacturing scale, media is shipped as a powder or granular (AGT) format, which is re-solubilised on site with sterile, pure water.
- At manufacturing scale no antibiotics (such as penicillin or streptomycin) are used in the medium.
 - For cost and allergy reasons, but does put the onus on all operators to ensure really good sterile techniques for cell culture.

Lecture Topics

A vertical diagram on the left side of the slide. It consists of four circles connected by a line. The top three circles are light gray, and the bottom circle is bright blue. Each circle is connected to a horizontal rectangular box on the right. The boxes are white with a dark blue border. The text inside the boxes is dark blue. The top box is labeled 'Cell Banking', the second 'Cell Culture Media', the third 'Defined Portion of Cell Culture Media', and the bottom, highlighted box is labeled 'Complex Portion of Cell Culture Media'.

Cell Banking

Cell Culture Media

Defined Portion of Cell Culture Media

Complex Portion of Cell Culture Media

Undefined Portion of Media

- Serum regarded as the black box of tissue culture media!
 - Usually makes up 5-20% of media – depending on cell line
- Serum contains many things including the following:
 1. **Hormones** – e.g. Growth hormone which increases growth, **Insulin** which promotes uptake of glucose
 2. **Growth factors** – e.g. Epidermal growth factor (EGF), nerve growth factor (NGF)
 3. **Proteins** can act as carriers for lipids and hormones e.g. α 2 macroglobulin inhibits trypsin
 4. **Lipids** – e.g. oleic acid, stearic and linoleic acid
 5. **Attachment factors** e.g. collagen, fibronectin and laminin
 6. **Trace elements** e.g. like Cu, Zn may act as enzyme co-factors and Se can inactivate free radicals

Advantages and Disadvantages of Serum

Advantages	Disadvantages
Allows the culture of many different cell types in non-defined media (Some cells will only grow in media with serum)	Quality of serum differs from batch to batch (Majority is foetal calf serum, also calf & horse serum)
Cells generally grow well and get high production	Supply may be restricted due to disease or drought
Serum offers protection against shear in bioreactor	Serum may act as a toxin to cells due to presence of viruses and mycoplasma
	Changing batches requires extensive testing
	Presents difficulties in downstream processing

Different Types of Media Available

1. Serum-free Media (SFM):

- These media do not require supplementation with serum but may contain discrete proteins or bulk protein fractions. In most cases the protein has been kept to a minimum and the media has been optimised for a specific cell type.

2. Protein-Free Media (PFM):

- These are a step closer to a **chemically-defined formulation** as compared to serum-free. However, they may still contain undefined components of animal or plant origin such as various hydrolysates that contribute low molecular weight peptides.

3. Chemically-Defined (CD) Media:

Contains no proteins, hydrolysates or components of unknown composition. All components have a known and well-defined chemical structure. This is the 'gold standard' of media available.

Different Types of Media Available

- **Generally cells can be adapted to grow in serum-free media**
 - Gradually over period of time, reduce amount of serum: can take several weeks or months
 - Can be more drastic and eliminate serum from media and select cells that grow in serum-free
- **Due to lack of serum require Pluronic®F68**
 - Pluronic is a non-ionic surfactant that is used to protect animal cells from damage caused by shear and the effects of sparging
 - Pluronic acts by creating an artificial layer around cells or is absorbed into membrane thus strengthening it
- Cells growing in serum-free media may not reach same cell densities as in serum containing media.

Advantages and Disadvantages of Serum-Free

Advantages	Disadvantages
Easier downstream processing : don't have to worry about getting rid of proteins in serum (e.g. albumins)	More expensive! : even though serum is costly, individual components of serum-free also costly
	Some cells more sensitive to shear in serum-free (use Pluronic F68)
Less problems with FDA as no risk of animal-derived products (viruses, prions)	Serum may act as a toxin to cells due to presence of viruses and mycoplasma
	Cells may not reach as high density as with serum
	Cells may be harder to freeze down

Protein-free Media

- **An increasing emphasis from regulatory authorities for the removal of animal-derived raw materials from antibody producing processes.**
- Number of groups have developed protein-free media:
 - Most of these are formulated for specific cell type : different for each cell type
- A number of successful formulations:
 - Protein-free medium for NS0 (hybridoma cell line) contains zinc instead of insulin, amino acids, ethanolamine, tropolone, nucleosides, glutamine
 - Protein-free medium for BHK : e.g. SMIF6
 - Many companies will use one media for growth of cells and then switch to another media for production

Protein-free Media

- Serum is a complex source of proteins, albumins, lipids, vitamins, growth factors, low molecular weight components, mineral ions etc.
- Serum proteins protect cells from shear, from protease attack & also facilitates transfer of vitamins through cells:
 - This reduces the risk of introducing adventitious agents and other contaminants. The main concerns relate to BSE (bovine spongiform encephalopathy) and CJD (Creutzfeldt -Jakob disease).
 - Makes process optimisation easier if poorly-defined additives such as serum and hydrolysates are avoided.
 - Improves process consistency
 - It reduces the protein load for downstream purification

Summary

- Cell culture media must contain all elements necessary to allow growth of cells including energy source, macro and micro nutrients
- Growth media usually contain 2 parts: defined part and non-defined portion:
 - Significant move towards fully defined media
- Defined part – key components:
 - Energy – glucose, pyruvate, Glu and Arg
 - Buffers – NaHCO_3 , CO_2
 - Amino acids – for building proteins, -source
 - Vitamins & essential fats
 - Inorganic ions – metals like Ca, Mg, Fe, Zn etc.
- Undefined portion - serum

Sample Questions

- Explain the structure of a cell bank system and the relationship between master and working banks.
- What steps are taken in the characterisation of cells for a cell bank?
- Discuss the following points: Characterisation tests include:
 - Viability (>90% viable) / Contamination free of microbes and a pure culture / Morphology / Isoenzyme profiles / Karyotype / DNA fingerprinting / Productivity - Titre
 - Master bank - failsafe - only access as a last resort. All cells stored derived from a single cell clone.
 - Use 1 master cell bank vial to expand and create multiple copies to act as working bank.

Sample Questions

- What are the major energy sources for mammalian cell contained in culture media? What concentration is each normally provided at? What impact would addition of additional sugar have on the medium – explain why?
- Discuss the following points:
 - Carbon sources: Glucose / Na-Pyruvate / Glutamine
 - Levels: Glucose – 4.5 to 6.0 g/L ; Na-Pyruvate – 0.11g/L; L-Gln – 0.584g/L; L-Arg – 0.084g/L
 - Excess glucose – upsets the osmolality of the medium
 - Low Glucose media = 324-333 mOsm
 - High Glucose media = 400 mOsm
 - If hypertonic, then cells risk dehydration and cell death.

Questions?

