



Reminder

 Have a look at the project titles and put your name down for your preferred project by 12th February





Learning Objectives

Discuss the complexity of biologics testing.

 Describe the different types of biologics testing; product characterisation, stability testing, in-process and product release testing.

List critical quality attributes of biologics

Identify key ICH guidelines for the quality control of biologics



Topics



Introduction to QC Testing

Key ICH Guidelines

What are Biopharmaceuticals?

- Biopharmaceuticals are medicinal products derived from biological sources
 - E.g. genetically engineered cells



 Bioprocessing is the manufacturing process for biopharmaceuticals






A World of Complexity

Traditional Medicines

Aspirin (180.2 Da) Atorvastatin (558.6 Da) Ibuprofen (206.3 Da)

Biopharmaceuticals






What makes proteins so complex?





What makes a protein?



 Proteins are primarily composed of molecular building blocks called 'amino acids'


 Many biopharmaceuticals also contain carbohydrates (sugars) called 'glycans'




What makes a protein?

 The amino acids are arranged in linear chains called 'polypeptides' by an enzyme complex in the cell

 Some proteins contain more than one polypeptide chain








What makes a protein?


- The amino acids are arranged in a specific sequence which is dictated by the **gene** for that protein
- The sequence is extremely important for the protein to take its biologically active shape





Post-translational Modification (PTM)

- Proteins can undergo 100's of additional modifications within the cell and also during bioprocessing
- These modifications usually significantly impact on stability and function
- Most important PTMs -Glycosylation




Pagel, Oliver, et al. "Current strategies and findings in clinically relevant post-translational modification-specific proteomics." *Expert review of proteomics* 12.3 (2015): 235-253.

Proteins are not static, rigid structures

The folding of the protein is not fixed!


 The protein will have a specific conformation ideal for biological interactions: 'active conformation'

- **Processing conditions** influence this folding
- Proteins can easily degrade if handled incorrectly





 A protein binding to its target is like a key fitting a lock. If the protein is the incorrect structure, it will not bind the target effectively.




 Preservation of this biologically active conformation is what drives the process design


Formulation and Stability Issues

- As traditional medicines degrade, major concern is loss of efficacy
- As biopharmaceuticals degrade, there can be loss of efficacy and potentially increased safety issues

- AggregationHydrolysis
- PrecipitationPhotolysis
- Fragmentation Deamidation (ASX,GLX)

- Oxidation (MET)
- Disulfide Scrambling
- Deglycosylation (Glycoproteins)








Process-related impurities must also be tackled!

Detecting and removing process related impurities represent

another challenge

- Host Cell Proteins
- Host Cell DNA
- Virus
- Residual Protein A






What about the process?






Chemical Synthesis: Follow the Recipe



Traditional pharmaceuticals are made through chemical synthesis




Many complex cell culture and purification procedures.





Upstream Variables



Upstream

Harvest

Downstream

Formulation and Fill



Cells


Mutations
Cell metabolism
Contaminants

Culture conditions

Medium composition Nutrient availability Bioreactor scale pH & temperature



Harvest Variables



Upstream

Harvest


Downstream

Formulation and Fill



Depending on harvest method used, could have increased levels of enzymes leaking from broken cells that will degrade the biopharmaceutical

Downstream Variables




Purification columns

Composition of purification solutions (buffers)

Exposure to pH extremes

Fill Finish Variables



Final formulation composition
Vials/syringe components
Contaminants from filling line equipment



Ensure Product Quality with Process Control

With so many variables involved, each step of the process must be controlled through **GMP activities** and monitored using **QC testing**





The Quality Management System (including QC testing) proves the consistency of the process

IRT ©



Topics



Introduction to QC Testing

Key ICH Guidelines



QC testing and Product Quality

Quality control testing aims to ensure the safety and efficacy of medicines

- Quality cannot be tested into the product; must be built into the process
- Thus, the process is monitored by testing the quality of:
 - All chemicals and solutions associated with the production of the medicine
 - The drug substance (DS)
 - The drug product (DP)





Quality Control testing in biopharma

Product characterisation

• Initial in-depth analysis of product. Set quality specifications (proven acceptable ranges)



Stability Testing

Determine the shelf life and storage conditions for product

- In-process testing

Testing bulk drug substance and associated materials

Lot release

 Ensure final product meets specifications as per product characterisation



Buffer/Media Prep

- ≽pH
- **≻**conductivity
- **≻**Osmolality
- **≻**Bioburden
- **≻**Endotoxin
- ➤ Water quality testing

Raw Materials

- ➤ Identity and purity
- **≻**HPLC
- ➤ Spectroscopy (FTIR, NMR)
- ➤ Residual moisture (KF Titration)
- **≻**TOC

Typical QC Testing in Bioprocessing

UPSTREAM

- ➤ Cell density and % viability
- ➤pH and gases
- ➤ Cell metabolites and substrates
- ➤ Biopharmaceutical titre (HPLC)
- ➤ Endotoxin and/or Bioburden samples (pre and post-innoculation, also in production)

HARVEST

- ➤ Total protein (protein assay)
- ➤ Biopharmaceutical titre (HPLC)
- ➤ Micro testing: Endotoxin and Bioburden (pre-filtration)
- External testing for mycoplasma and virus pre-filtration

DOWNSTREAM

- ➤ Total protein by UV after UF
- ▶pH/conductivity after DF
- ➤ Biopharmaceutical titre after every chromo step
- >Stability indicating assays (HPLC, electrophoresis)
- ➤ Bioburden before each step, i.e. pre filtration, pre chromo etc.
- > Endotoxin performed after each step


Formulated DS & DP

- ➤ Purity, integrity, identity (HPLC, electrophoresis)
- ➤ Total protein by UV or Bradford/BCA
- ➤ Potency (cell-based assay/bioassay)
- ➤ Process related impurities (Western blot, ELISA, HPLC, SPR)
- ➤ Product related impurities (HPLC, electrophoresis)
- ➤ Peptide mapping (HPLC)
- ➤ Glycan analysis (HPLC/Capillary electrophoresis)
- ➤ Charge variants (IEF, HPLC)
- ➤ Stability testing
- ≻pH
- **≻**conductivity
- **≻**Osmolality
- **≻**Appearance
- ➤ Leachables and extractables
- **≻**Bioburden
- **≻**Endotoxin
- ➤ Extractable volume




Critical Quality Attributes (CQAs)

- 'A CQA is a physical, chemical, biological, or microbiological property or characteristic that should be within an appropriate limit, range, or distribution to ensure the desired product quality.' ICH Q8 (R2)
- E.g. molecular weight, potency, charge variants, glycoforms, process and product related impurities
- Any property that affects:
 - Biological Activity
 - PK/PD
 - Immunogenicity
 - Safety





Quality Control testing ensures all CQAs are within specification



A specification is a **proven acceptable range** (PAR) for a given QC test result

Ramanan, Sundar, and Gustavo Grampp. "Drift, Evolution, and Divergence in Biologics and Biosimilars Manufacturing." BioDrugs (2014): 1-10.



How does the manufacturer know what tests to use?

- Manufacturer is free to use any scientifically valid test for testing of product in their license application
- Information on generally acceptable types of testing:
 - Pharmacopoeias (EP, USP, etc.)
 - FDA Guidance documents
 - ICH Guidance documents
 - Scientific literature, etc.
- Also possible to use alternative tests to those prescribed by biological product standards: 21 CFR 610.9/ICH, if justified
- Once approved, must continue to use the approved tests for all subsequent batches



Topics



Introduction to QC Testing

Key ICH Guidelines

ICH




 International Council for Harmonisation of Technical Requirements for Registration of Pharmaceuticals for Human Use

- Unique harmonisation project involving the regulators and industry of:
 - US, EU and Japan
 - WHO, Canada, and EFTA are observers

ICH



ICH Provides 60 Guidelines on technical requirements under 4 Working Groups



http://www.ich.org/products/guidelines/quality/article/quality-guidelines.html

ICH Guidelines for Quality of Biopharmaceuticals



Quality Guidelines


Q5A - Q5E Quality of Biotechnological Products				
Code	Document Title	Previously code		
Q5A(R1)	Viral Safety Evaluation of Biotechnology Products Derived from Cell Lines of Human or Animal Origin	Q5A		
Q5B	Analysis of the Expression Construct in Cells Used for Production of r-DNA Derived Protein Products			
Q5C	Stability Testing of Biotechnological/Biological Products			
Q5D	Derivation and Characterisation of Cell Substrates Used for Production of Biotechnological/Biological Products			
Q5E	Comparability of Biotechnological/Biological Products Subject to Changes in their Manufacturing Process			






ICH Topic Q 6 B

Specifications: Test Procedures and Acceptance Criteria for Biotechnological/Biological Products




Product-related impurities




Disulfide bridges




Intact mass
Higher Order Structure
Potency
Binding properties



Glycan profile









Vial inspection






















Immunoassays





Microbial testing







Quality Guidelines


Q5A - Q5E Q	uality of Biotechnological Products	
Code	Document Title	Previously coded
• Q5A(R1)	Viral Safety Evaluation of Biotechnology Products Derived from Cell Lines of Human or Animal Origin	Q5A
▶ Q5B	Analysis of the Expression Construct in Cells Used for Production of r-DNA Derived Protein Products	
Q5C	Stability Testing of Biotechnological/Biological Products	
Q5D	Derivation and Characterisation of Cell Substrates Used for Production of Biotechnological/Biological Products	
▶ Q5E	Comparability of Biotechnological/Biological Products Subject to Changes in their Manufacturing Process	





Stability Testing

Key part of QC Program and conducted throughout the product lifecycle



To understand how the **quality** of a DS/DP **changes over time** under the influence of **different environmental conditions**, and to determine degradation pathways

Provides information with regard to DS/DP storage conditions, final formulation/packaging and shelf life. This information is necessary for regulatory approval

Ensure continuous safety and efficacy for the patient

SGS, "A guide to biologic stability testing" Accessed April 2002 https://www.sgs.com/en/news/2015/10/a-guide-to-biologic-stability-testing

Blessy, M., et al. (2014) "Development of a forced degradation and stability indicating studies of drugs - a review" Journal of Pharmaceutical Analysis 4(3) pp 159-165

Product must be completely characterised, and CQAs determined before stability studies begin



Stability Testing: What conditions do we look at?

Temperature

 Storage temperature and above

Humidity

 Can sometimes be omitted where humidity protecting containers are used (Climate zones I and II)

Light

Studies to demonstrate the effects of light
Case by base basis

Container/ Closure Systems

- Inverted, horizontal and upright studies
- Effects of closure

Reconstitution

- Stability after reconstitution
- Storage conditions

The purpose of stability testing is to:

provide evidence on how the quality of a drug substance or product varies over time, under the influence of a variety of environmental factors

ICH Q5C: Stability Testing of Biotechnological/Biological Products ICH Q1A(R2): Stability Testing of New Drug Substances and Products ICH Q1B: Photostability Testing of New Active Substances and Medicinal Products





Quality Guidelines

Q5A - Q5E Q	uality of Biotechnological Products	
Code	Document Title	Previously coded
Q5A(R1)	Viral Safety Evaluation of Biotechnology Products Derived from Cell Lines of Human or Animal Origin	Q5A
• Q5B	Analysis of the Expression Construct in Cells Used for Production of r-DNA Derived Protein Products	
• Q5C	Stability Testing of Biotechnological/Biological Products	
Q5D	Derivation and Characterisation of Cell Substrates Used for Production of Biotechnological/Biological Products	
▶ Q5E	Comparability of Biotechnological/Biological Products Subject to Changes in their Manufacturing Process	





Changes in CQAs: Drift and Evolution in Biologics

Changes in the biologic can be a result of:

 Drift: a result of unknown deviations in the manufacturing process

 Evolution: known changes in the manufacturing process (equipment change, scaling, raw materials suppliers, etc.)








Photo: https://www.wuxibiologics.com/scale-out-vs-scale-up-biomanufacturing/

Ramanan, Sundar, and Gustavo Grampp. "Drift, Evolution, and Divergence in Biologics and Biosimilars Manufacturing." *BioDrugs* (2014): 1-10.

4



The goal of ICH Q5E

- The goal of this comparability exercise is to ensure the quality, safety and efficacy of a drug product produced by a changed manufacturing process
- Demonstration of comparability does not necessarily mean that CQAs of pre-change and post-change product are identical
 - They are highly similar and that existing knowledge is sufficiently predictive to ensure that any differences in CQAs have no adverse impact



Gather important information



Identify impact of changes on product quality



Select relevant analytical methods



Determine acceptance criteria

- List of CQAs
- Description of product change and rationale for change
- Historical data and product characterisation data

- What might be effects on CQAs?
- What might be effects on inprocess controls?

 Which analytical methods are most relevant to detect potential change in corresponding CQA Based on side by side comparison or statistical comparison




Example of a Comparability Plan

Process Change	Affected CQA	Impact Assessment	Analytical Method	Acceptance Criteria
	Residual HCP	Scale-up is expected to produce more biomass hence more residual HCP	HCP ELISA	≤ 50ppm
Scale-up of cell	Glycosylation Profile	Scale up might affect cell	Oligosaccharide mapping	Side by side profile
culture Scale-up 10,000L	diycosylation Frome	line growth parameters possibly leading to a change in the glycosylation pattern	Oligosaccharide mapping	comparison
	Isoform profile	Scale up might affect cell	HPLC	Main peak ≥ 90%
		line growth parameters possibly generating a different isoform profile		Pre-peak < 3% Post-peak < 2%



Change Risk and Data Requirements



The nature of a manufacturing change determines the amount and type of supporting data required to evaluate comparability



Topics



Introduction to QC Testing

Key ICH Guidelines



Thank You

